


BAUERFEIND.FOR A BETTER QUALITY OF LIFE





Bauerfeind AG is one of the world's leading manufacturers of medical aids. With our products and health care concepts, we enable people to get moving and thus play an important role in helping them stay healthy.

We have been developing compression products for phlebological and lymphatic conditions for decades, approaching our production and service with untiring dedication and our development with ingenuity.

Thanks to our custom-made VenoTrain curaflow flat knit products, we can offer our lymphatic patients perfectly tailored solutions, with six different colors to choose from (cream, caramel, anthracite, black, navy, and bordeaux) and more than 1,000 different ways of combining products and additional features. The beneficial therapeutic features of the range include the knit's high microfiber content, the flat sensitive stitching, and special comfort zones such as the extra-soft area behind the knee.

The flat-knit compression stocking VenoTrain delight is indicated for use with advanced venous disorders and 1st stage lymphedema. The breathable double-stretch knit is latexfree, skin-friendly, and easy to put on and take off. The reduced compression of the gripRelief structure on all supportive bands prevents unpleasant constriction.

Our circular knit VenoTrain impuls compression stocking combines a wide range of benefits: Its breathable microfiber and firm compression thread offers unique wearing comfort with a strong effect. Only very little movement is needed for the stocking to exert the required therapeutic pressure.

CONTENTS

2 Bauerfeind lymph expertise



LEG & FOOT

7 Care system

VenoTrain impuls

9 Standard & custom-made product VenoTrain delight



- 11 Custom-made product
 - + additional features

VenoTrain curaflow

- 15 Custom-made product
 - + additional features



ARM & HAND

- 23 Care system
 - VenoTrain soft
- 25 Standard product
 - VenoTrain curaflow
- 27 Custom-made product + additional features



THORAX VenoTrain curaflow

- 33 Custom-made product + additional features
- 36 ACCESSORIES

SERVICE & ORDERING INFORMATION

- 38 Alterations | Repairs
- 40 Sizing systems
- 41 Order pads

A RELIABLE COMBINATION ACCORDING TO USERS*



COMFORTABLE!

THE MATERIAL

90 percent of users rated the VenoTrain curaflow material as good or excellent.



WELL-FITTING!

THE FIT

91 percent of those surveyed felt that their VenoTrain curaflow compression stockings fitted well or very well.



CONVENIENT!

PUTTING ON THE PRODUCT

82 percent of users also found their VenoTrain curaflow products easy or very easy to put on.

^{*} Source: Results from a patient survey assessing treatment with VenoTrain curaflow compression products, Bauerfeind AG – March to August 2014





DESIGNED FOR EVERYDAY LIFE

- Smooth microfibers make the stockings convenient to put on and comfortable to wear
- Constant pressure and an ideal pressure profile with a tangible effect throughout the time worn, pleasantly soft, fine double-stretch knit
- Perfect fit for individual body dimensions thanks to anatomically tailored stitched products
- Breathable and moisture-regulating knit for excellent wearing comfort

VenoTrain curaflow ADDITIONAL FEATURES

- Extra-soft behind the knee for problem-free mobility
- Intelligent longitudinal structure for optimal lymph transportation

UNCOMPROMISING PRODUCT QUALITY





Our VenoTrain compression stockings cover all areas of application associated with modern compression therapy. They offer impressive, clinically proven effectiveness and unique material properties.

Thanks to the wide range of colors available, the compression products blend seamlessly into their users' everyday lives and can be combined with various fashion styles.

As a preventive measure prior to treatment of lymphedema: The solutions, which are tailored to specific indications, include a wide variety of circular and flat knit styles that are available as custom-made or standard products.

All VenoTrain products conform to the requirements of the German Medical Devices Act (MPG) and bear the RAL quality mark. As a seal of quality, the Standard 100 by Oeko-Tex® guarantees testing for harmful substances.







COMPREHENSIVE LYMPH EXPERTISE BAUERFEIND ACADEMY

It is important to regularly train the correct measurement of flat knit products in order to ensure the best possible care for patients.

In cooperation with your local representative, Bauerfeind Academy offers individualized training sessions on lymphology.

Please contact your local organisation for further details.







LYMPHOLOGY TRAINING CONTENTS:

- Medical basics of lymphology and compression therapy
- Complex physical decongestive therapy
- · Indication-based treatment with Bauerfeind flat knit products
- Practical measurement and fitting of flat knit products on the leg, arm, hand and thorax including additional features such as toe caps, pads, etc.
- Practical work involving direct patient contact (measuring and fitting)



CARE SYSTEM LEG & FOOT

STAGES OF LYMPHEDEMA & LIPEDEMA

Stage 0 **LATENCY STAGE**

The lymphatic system is damaged, but there are still no visible or palpable edemas.



Stage 1
SPONTANEOUSLY
REVERSIBLE STAGE

The tissue is soft and can be pressed down gently with a visible indentation. The lymphedema subsides when the leg or arm is elevated.

The surface of the skin is smooth, the subcutaneous adipose tissue is thickened, and the fat structure is finely knotted.



Stage 2
SPONTANEOUSLY
IRREVERSIBLE STAGE

The tissue is hardened and it is either difficult or impossible to press down on it and create an indentation. There may be pain and restriction of movement.

The surface of the skin is uneven and the fat structure is coarsely knotted.



Stage 3 **ELEPHANTIASIS STAGE**

A significant increase in edema size, sometimes with extreme dimensions.

The tissue is also tougher and harder. There are also lobes of fat causing extensive deformation

PRODUCT RECOMMENDATION

VenoTrain impuls stages 0-1

VenoTrain delight stages 0-2

VenoTrain curaflow stages 0-3

VenoTrain impuls COMPRESSION CCL 2-3

- Soft microfiber, smooth and breathable
- · Short-stretch elastic
- · Circular knit

VenoTrain delight COMPRESSION CCL 2-4

- Latex-free and gentle on the skin
- · Hard-wearing
- · Extra-flat seams, flat knit

VenoTrain curaflow COMPRESSION CCL 1-4

- Soft microfiber, smooth and breathable
- Extra-flat seams, flat knit

COMPRESSION CLASSES

CCL 1: 18-21 MMHG

- · Light compressio
- · Light surface effect
- Treatment of childhood lymphedema, edema caused by inactivity

CCL 2: 23-32 MMHG

- · Moderate compression
- Moderate surface effect
- · Light deep-vein effect
- Treatment of primary and secondary lymphedema of the arm and leg (stages 1 and 2), lipedema, lipolymphedema

CCL 3: 34-46 MMHG

- . Firm compression
- Surface and deep-vein effect
- Treatment of primary and secondary lymphedema of the arm and leg (stages 2 and 3), lipedema, lipolymphedema

CCL 4: FROM 49 MMHG

- · Extra-firm compression
- · Strong deep-vein effect
- Treatment of primary and secondary lymphedema of the leg (stages 2 and 3), lipedema, lipolymphedema





VenoTrain® impulsFOR ADVANCED VEIN PROBLEMS

Produces a particularly deep action by exerting high working pressure with even the slightest muscle movements. Easy to put on and breathable thanks to microfibers.

Highly effective for advanced vein problems, post-thrombotic syndrome, follow-up care for healed ulcers, and therapy in the early stages of lymphedema.

COLORS









FEATURES

- The short-stretch knit produces a high working pressure with even the slightest muscle movements
- Easy to put on thanks to microfiber material
- Skin-friendly microfiber material makes the product breathable and ensures climate regulation
- Highly differentiated V6+ (I–VI) sizing system

DETAILS

SIZING SYSTEM V6+ (I–VI), custom-made MATERIAL CONTENT

- · Ccl 2: 80 % polyamide, 20 % elastane
- · Ccl 3: 75 % polyamide, 25 % elastane

COMPRESSION CLASSES

- · Ccl 2 (23–32 mmHg)
- · Ccl 3 (34-46 mmHg)

TOE open or closed

CARE INSTRUCTIONS



STANDARD & CUSTOM-MADE SOLUTIONS

TOP BAND EDGING AD AG AGHB AT ATH ATX Top band Silicone dots Garter belt



VenoTrain® delight FOR INTENSIVE THERAPY

Available as a flat-knit leg stocking with skinfriendly, latex-free materials and exceptional fit.

For the treatment of advanced venous and lymphatic symptoms. The VenoTrain delight is available in a wide variety of designs.

COLORS









FEATURES

- · Deep effect on the venous system
- · Breathable double-stretch fabric
- · Special gripRelief® structure at the top ensures the stocking remains comfortably and securely in place
- · Extra-flat seams

Viscoelastic

SIZING SYSTEM custom-made **MATERIAL CONTENT**

- · Ccl 2: 80 % polyamide, 20 % elastane
- · Ccl 3: 79 % polyamide, 21 % elastane
- · Ccl 4: 75 % polyamide, 25 % elastane

COMPRESSION CLASSES

- · Ccl 2 (23–32 mmHg)
- · Ccl 3 (34-46 mmHg)
- · Ccl 4 (from 49 mmHg)

TOE open, diagonal cut or closed

CARE INSTRUCTIONS





CUSTOM-MADE SOLUTIONS

TOP BAND EDGING	AD	AF	AG	AGHB	AT	BD	BG	CG	CT	DG	ET
Top band	•				•	•	•	•	•	•	•
Silicone dots	•	•	•				•	•		•	
Garter belt				•							
ADDIT. FEATURE											
T-heel											

VenoTrain® delightADDITIONAL FEATURES FOR LEG & FOOT





01 FOOT WITH DIAGONAL CUT PATIENT'S SITUATION

- · Edematous dorsum of the foot
- · Toe misalignments
- · Skin and nail conditions in the toe area

EFFECT

- Required pressure profile at the metatarsophalangeal joints
- · Anatomical fit

FOOT STRAIGHT

PATIENT'S SITUATION

- Edematous dorsum of the foot
- Toe and foot misalignments
- Narrow foot tapering toward the heel

EFFECT

- Prevention of foot and/or heel slipping backward
- Very well tolerated soft toe

03 T-HEEL

PATIENT'S SITUATION

- · Skin folds at Y
- Pronounced anatomical heel

EFFECT

- No creasing at the heel or instep
- · Improved therapeutic gradient in the ankle region
- · Better joint functionality







0.4

04 SILICONE DOTTED BAND

PATIENT'S SITUATION

 Generally required on all compression stockings for fastening in place

EFFECT

- Holds the compression stocking in place
- · User-friendly

05 PAD

PATIENT'S SITUATION

· Edema and hardening in the ankle region

EFFECT

- Increase in pressure at specific points
- · Prevention of edemas and hardening




0.5













VenoTrain® curaflow FOR TREATING LYMPHEDEMA AND LIPEDEMA

The high microfiber content and the innovative way these are processed make this item unique. Extra-flat sensitive stitching and special comfort zones ensure this flat knit product provides excellent wearing comfort.

COLORS



FEATURES









DETAILS

- Smooth microfibers make the product convenient to put on
- Skin-friendly microfiber material makes the product breathable and ensures climate regulation
- · Extra-soft toe
- Anatomical, extra-soft area behind the knee for unrestricted freedom of movement
- Custom-made solutions to suit individual requirements

SIZING SYSTEM custom-made

MATERIAL CONTENT

- · Ccl 1–2: 73 % polyamide, 27 % elastane
- · Ccl 3: 70 % polyamide, 30 % elastane
- \cdot Ccl 4: 65 % polyamide, 35 % elastane

COMPRESSION CLASSES

· Ccl 1 (18–21 mmHg) to Ccl 4 (from 49 mmHg) **TOE** open, diagonal cut,

closed, closed at an angle

CARE INSTRUCTIONS



CUSTOM-MADE SOLUTIONS



VenoTrain® curaflow ADDITIONAL FEATURES FOR THE LEG & FOOT

Thanks to the custom-made solutions available. VenoTrain curaflow products can be individually adapted to suit each wearer's requirements. The product features shown here are the most popular additional options. Our extensive range of custom-made solutions offers a total of over 1,000 possible product combinations. Please refer to our current price list for more additional features





FOOT CLOSED AT AN ANGLE **PATIENT'S SITUATION**

· Edematous dorsum of

- the foot
- · Toe and foot misalignments
- Narrow foot tapering toward the heel

EFFECT

- · Required pressure profile at the metatarsophalangeal
- Anatomical fit
- Prevention of foot/heel slipping backward
- Very well tolerated soft toe

FOOT STRAIGHT AND OPEN PATIENT'S SITUATION

- · Edematous dorsum of the foot with a relatively straight toe line
- Foot not tapered

EFFECT

- Prevention of foot/heel slipping backward
- Very well tolerated soft toe

FOOT WITH DIAGONAL CUT PATIENT'S SITUATION

- · Edematous dorsum of the
- Toe misalignments
- Skin and nail conditions in the toe area

EFFECT

- Required pressure profile at the metatarsophalangeal ioints
- · Anatomical fit

90° HEEL **PATIENT'S SITUATION**

· Skin folds at Y

- · Pronounced anatomical heel
- · Wheelchair user

- · No creasing at the heel or instep
- · Improved therapeutic gradient in the ankle region
- · Better joint functionality

INSTEP WIDENING PATIENT'S SITUATION

- · Bulges at Y
- · Lymphedema (end of stage 2/beginning of stage 3)

EFFECT

- Reduced pressure on the dorsiflexor tendon
- Improved therapeutic gradient on the dorsum of the foot

HALLUX

PATIENT'S SITUATION

· Bunion (hallux valgus)

EFFECT

- · Optimum fit in the region of the big toe
- · Pressure relief















07 KNEE ZONE (STANDARD FEATURE ON ALL LEG PRODUCTS WITH A KNEE) PATIENT'S SITUATION

Generally for lymph and lipedema patients

EFFECT

- Extra-soft behind the knee to reduce creasing
- Anatomical fit for better joint functionality

08 SLANTED EDGING PATIENT'S SITUATION

- Generally for lymph and lipedema patients
- · Tapered anatomical shape
- · Soft skin tissue

EFFECT

- Improved therapeutic gradient at the ends of the stockings (AD and AG possible)
- · Anatomical fit

09 SILICONE DOTTED BAND PATIENT'S SITUATION

 Generally required on all compression stockings for fastening in place

EFFECT

Holds the compression stocking in place

STRIPED SILICONE TOP BAND PATIENT'S SITUATION

 Offers a stronger hold; suitable for firmer edemas and large leg circumferences

EFFECT

Holds the compression stocking in place

11 SILICONE TOP BAND SECTION

PATIENT'S SITUATION

- Pronounced tapered leg shape
- · Very soft tissue

EFFECT

· Extra grip









Sino



12 **GARTER BELT PATIENT'S SITUATION**

· Cannot tolerate silicone top band

EFFECT

· Holds the compression stocking in place

ZIP FASTENING

PATIENT'S SITUATION

· Secondary conditions such as rheumatism, arthritis, incontinence, etc.

EFFECT

· Makes it easier to put on and take off the product so that it fits correctly

ATU GUSSET PATIENT'S SITUATION

· Pregnancy

EFFECT

- · Improved fit in the panty
- · Reduced compression in the abdominal region

15 **STANDARD MEN'S GUSSET**

PATIENT'S SITUATION

· Pantyhose for men









16 LOOPED/PULL-THROUGH **ELASTIC BAND (STANDARD)**

· Fluctuating waist circumference

PATIENT'S SITUATION

- · Pregnancy
- · Asthmatic

EFFECT

- · Adjustable width in the waist area
- · Adapts to changing body shape

17

WAIST BAND

PATIENT'S SITUATION

· For all pantyhose products

EFFECT

· Perfect fit without being bulky

INTEGRATED PAD PATIENT'S SITUATION

· Edema and hardening in the ankle region

EFFECT

- · Increase in pressure at specific points
- · Prevention of edemas and hardening

INTEGRATED LYMPH PAD PATIENT'S SITUATION

- · Locally aggravated edema or hardening
- · Variations from standard anatomical shape

EFFECT

- · Increase in pressure at specific points
- · Prevention of edemas and hardening

20 **TOE CAP**

PATIENT'S SITUATION

- · Forefoot and toe lymphedema
- · Edema of the dorsum of the

EFFECT

· Prevention of edema in the forefoot and toe area

SEWN-ON TOE LOOPS PATIENT'S SITUATION

- · Forefoot and toe lvmphedema
- · Edema of the dorsum of the

EFFECT

- · Prevention of edema in the forefoot and toe area
- · No overlapping of the stocking and toe section

















CARE SYSTEM ARM & HAND

STAGES OF LYMPHEDEMA & LIPEDEMA

Stage 0 **LATENCY STAGE**

The lymphatic system is damaged, but there are still no visible or palpable edemas.

Stage 1

SPONTANEOUSLY REVERSIBLE STAGE

The tissue is soft and can be pressed down gently with a visible indentation.

The surface of the skin is smooth, the subcutaneous adipose tissue is thickened, and the fat structure is finely knotted.

Stage 2

SPONTANEOUSLY IRREVERSIBLE STAGE

The tissue is hardened and it is either difficult or impossible to press down on it and create an indentation.

The surface of the skin is uneven and the fat structure is coarsely knotted.

Stage 3 **ELEPHANTIASIS STAGE**

A significant increase in edema size, sometimes with extreme dimensions.

The tissue is also tougher and harder. There are also lobes of fat causing extensive deformation.

PRODUCT RECOMMENDATION

VenoTrain soft stages 0-1

VenoTrain curaflow stages 0-3

VenoTrain soft Arm COMPRESSION CCL 1, CCL 2

- Comfort elbow for greater freedom of movement, wearing comfort, and a reduced pressure sensation
- Noticeable effect throughout the entire wearing period
- Gentle on the skin, soft, and breathable

VenoTrain curaflow COMPRESSION CCL 1-3

- Soft microfiber, smooth and breathable
- Extra-flat seams, flat knit

COMPRESSION CLASSES

CCL 1: 18-21 MMHG

- . Light compression
- · Light surface effect
- Treatment of childhood lymphedema, edema caused by inactivity

CCL 2: 23-32 MMHG

- Moderate compression
- Moderate surface effect
- · Light deep-vein effect
- Treatment of primary and secondary lymphedema of the arm and leg (stages 1 and 2), lipedema, lipolymphedema

CCL 3: 34-46 MMHG

- . Firm compression
- · Surface and deep-vein effect
- Treatment of primary and secondary lymphedema of the arm (stages 2 and 3), lipedema, lipolymphedema



VenoTrain® soft EVERYDAY MOBILITY

Seamless with an extra-soft inner elbow area. This facilitates everyday mobility with complete medical efficacy.

COLORS







DETA

- Comfort elbow for greater freedom of movement, wearing comfort, and
- · a reduced pressure sensation
- Noticeable effect throughout the entire wearing period
- · Gentle on the skin, soft, and breathable
- **SIZING SYSTEM** VenoTrain soft Arm (S–XL) **MATERIAL CONTENT**
- · Ccl 1: 74 % polyamide, 26 % elastane
- · Ccl 2: 69 % polyamide, 31 % elastane

COMPRESSION CLASSES

- · Ccl 1 (18–21 mmHg)
- · Ccl 2 (23–32 mmHg)

CARE INSTRUCTIONS



STANDARD

FEATURES

EDGING/ TOP BAND	CG
Top band	•
Silicone dots	•













VenoTrain® curaflow FOR TREATING LYMPHEDEMA AND LIPEDEMA **ARM & HAND**

The high microfiber content and the innovative way these are processed make this item unique. Extra-flat sensitive stitching and special comfort zones ensure this flat knit product provides excellent wearing comfort.

COLORS



FEATURES











- Smooth microfibers make the product convenient to put on
- Skin-friendly microfiber material makes the product breathable and ensures climate regulation
- · Anatomically shaped, extra-soft inner elbow area for greater freedom of movement, wearing comfort, and a reduced pressure sensation
- · Custom-made solutions to suit individual requirements

SIZING SYSTEM custom-made **MATERIAL CONTENT**

- · Ccl 1: 77 % polyamide, 23 % elastane
- · Ccl 2: 77 % polyamide, 23 % elastane
- · Ccl 3: 73 % polyamide, 27 % elastane Seamless hand
- · Ccl 1: 87 % polyamide, 13 % elastane
- · Ccl 2: 85 % polyamide, 15 % elastane
- · Ccl 3: 86 % polyamide, 14 % elastane Hand with seam
- · Ccl 1: 73 % polyamide, 27 % elastane
- · Ccl 2: 73 % polyamide, 27 % elastane
- · Ccl 3: 70 % polyamide, 30 % elastane

COMPRESSION CLASSES

· Ccl 1 (18–21 mmHg) to Ccl 3 (34–46 mmHg)

CARE INSTRUCTIONS



CUSTOM-MADE SOLUTION

















VenoTrain® curaflowADDITIONAL FEATURES FOR THE ARM & HAND

Thanks to the custom-made solutions available, VenoTrain curaflow products can be individually adapted to suit each wearer's requirements. The product features shown here are the most popular additional options. Our extensive range of custom-made solutions offers a total of over 1,000 possible product combinations. Please refer to our current price list for more additional features.





OT ARMZONE (STANDARD FOR ALL ARM SLEEVES) PATIENT'S SITUATION

Generally for lymph and lipedema patients

EFFECT

- Comfort elbow to reduce creasing
- Anatomical fit for better joint functionality

02 SLANTED EDGING PATIENT'S SITUATION

- Generally for lymph and lipedema patients
- · Tapered anatomical shape
- · Soft skin tissue

EFFECT

- Improved therapeutic gradient at the ends of the product
- · Anatomical fit
- Prevents circular constriction





03 SILICONE DOTTED BAND PATIENT'S SITUATION

 Required for fastening in place where the tissue is soft

EFFECT

- Holds the compression stocking in place
- · User-friendly

04 SILICONE TOP BAND SECTION

PATIENT'S SITUATION

- · Pronounced tapered leg shape
- · Very soft tissue

EFFECT

- · Extra grip
- · Prevents slipping

STRIPED SILICONE TOP BAND

PATIENT'S SITUATION

· Required for fastening in place where the tissue is soft

EFFECT

 Holds the compression stocking in place









06 ZIP FASTENING PATIENT'S SITUATION

 Secondary conditions such as rheumatism, arthritis, incontinence, etc.

EFFECT

 Makes it easier to put on and take off the product so that it fits correctly

07 INTEGRATED LYMPH PAD PATIENT'S SITUATION

- Locally aggravated edema or hardening
- · Variations from standard anatomical shape

EFFECT

- Increase in pressure at specific points
- · Prevention of edemas and hardening

08 SEAMLESS GLOVE PATIENT'S SITUATION

- Mild to moderate edema of the finger
- Mild to moderate edema on the back of the hand

EFFECT

 Prevention of edemas in fingers and on the back of the hand

09 GLOVE WITH SEAM PATIENT'S SITUATION


- Moderate to pronounced edema of the finger
- Moderate to pronounced edema on the back of the hand

EFFECT

 Prevention of edemas in fingers and on the back of the hand















VenoTrain® curaflow FOR TREATING LYMPHEDEMA AND LIPEDEMA **THORAX**

Lymph relief bras are generally insufficient for treating breast and chest wall edemas and individual solutions are often required: In these cases, thorax compression therapy with the VenoTrain curaflow Body can be the key to a better quality of life.

The seam-free chest area fits snugly and comfortably against the chest and exerts even compression on the lymphedema without causing any constriction or rubbing.

COLORS



ADDITIONAL FEATURES

- · Without arm sleeve
- · Arm sleeve with compression
- · Arm sleeve without compression
- · Central zip fastening at the front
- · Hook fastening
- · Sewn-on arm sleeve

FEATURES







DETAILS

· Smooth microfibers make the product convenient to put on

- Skin-friendly microfiber material makes the product breathable and ensures climate regulation
- Custom-made solutions to suit individual requirements

SIZING SYSTEM custom-made **MATERIAL CONTENT**

· Ccl 1: 73 % polyamide, 27 % elastane

COMPRESSION CLASSES

· Ccl 1 (18–21 mmHg)

CARE INSTRUCTIONS



CUSTOM-MADE SOLUTIONS











VenoTrain® curaflowADDITIONAL FEATURES FOR THE THORAX



01 SEAMLESS CHEST AREA PATIENT'S SITUATION

Standard feature on all thorax products

FFFFCT

- Prevents pressure points and chafe marks in the chest area
- · Supports the chest

02 ARM SLEEVE WITH OR WITHOUT COMPRESSION

PATIENT'S SITUATION

· Edema in the chest and arm area

EFFECT

 Prevention of edemas in the shoulder and upper arm area

03 EPITHESIS POCKET PATIENT'S SITUATION

· Use of breast prostheses

EFFECT

 Optimum grip when using breast prostheses

ZIP FASTENING PATIENT'S SITUATION


 Secondary conditions such as rheumatism, arthritis, incontinence, etc.

EFFECT

 Makes it easier to put on and take off the product so that it fits correctly

01









VenoTrain® ACCESSORIES

THE PERFECT ADDITION FOR EVERYDAY USE

The handy extras from the VenoTrain accessories range are useful for ensuring that patients feel comfortable using and wearing their VenoTrain compression stockings. They provide effective support for treatment, make it easier to handle the products, and care for the compression stockings as well as your customers' skin.

05 VenoTrain® lotion

sensitive skin

protection

developed in collaboration with

06 Bauerfeind care foam

07 Special detergent

Sebapharma; contains chamomile extract and is specially designed for people with

contains panthenol to protect and soothe particularly stressed and sensitive skin

cleans gently, ensuring careful treatment of your products, and providing lasting

01 Rubber gloves

help with putting the compression stockings on correctly

02 Measurement bag

for assistance when carrying out measurements

03 Starter kit

including gloves, a foot file, lingerie bag, non-skid soles, an exercise band, and a loyalty card

04 Washing bag

to protect the compression stockings during every machine wash

















VenoTrain® glider PRACTICAL DONNING AID

The VenoTrain glider is a practical aid for putting on and taking off medical compression stockings. It can be prescribed by a physician.

The lightweight, hard-wearing material makes the stocking glide effortlessly over the arm or leg. It is suitable for compression stockings with an open or closed toe.

Saving space and time, it is the ideal companion when you're out and about.

FEATURES

- · Extremely lightweight, robust, and tear-resistant
- · Easy to handle
- · Ideal for travel













ALTERATIONS | REPAIRS ADDITIONAL SERVICES INDIVIDUAL CUSTOMER SERVICE

ALTERATIONS | REPAIRS

STANDARD

		soft Arm	impuls Leg
Top of	Replace silicone top band	•	•
stocking	Top band to silicone top band	•	•
	Silicone top band	•	•
	Remove hip section and sew on silicone top band		•
AGHB	Replace hip strap		•
Shorten stocking	AD/CG by up to 3 cm	•	•
	AG by up to 5 cm		•
Repair,	Hole	•	•
as far as possible	Run	•	•
	Seam	•	
	Loose thread	•	•

REORDERING

CUSTOM-MADE SOLUTION

	Unit	impuls	delight
Knitted toe area Lengthen or shorten	Item	•	
Toe Seam on the outside	Item	•	•
Toe Make wider, narrower, anatomical	Item	•	•
Panty area Front/back – two seams	Item	•	
AT with men's gusset (without fly or with longitudinal slit)	Item	•	
AT with extra-wide "saddlebag" gusset	Item	•	
Sew in a foam pad	Item		•

SERVICE & ORDERING INFORMATION

ALTERATIONS | REPAIRS CUSTOM-MADE SOLUTION

ALTERATIONS | REPAIRS VenoTrain® curaflow


		impuls Leg	deligh Leg
Top of stocking	Replace silicone top band	•	•
	Top band to silicone top band	•	•
	Silicone top band narrower or wider	•	•
	Remove hip section and sew on silicone top band	•	•
AT	Shorten pantyhose top	•	•
	Lengthen pantyhose top (excluding ATX)	•	
	Take in gusset (excluding ATX)	•	
	Widen gusset (excluding ATX)	•	
	Extra-wide "saddlebag" gusset	•	
	Alter AT to AG with silicone top band	•	
AGHB	Replace hip strap/ shoulder strap	•	•
Shorten	AD by up to 3 cm	•	•
stocking	AF by up to 4 cm		•
	AG by up to 5 cm	•	•
Redye from	ı light to dark*	•	•
Pads**	Sew in customer's own pad		•
	Sew in a foam pad		•
Toe alteration	Alter closed toe to open toe	•	•
	Shorten open toe	•	•
Repair,	Hole	•	•
as far as	Run	•	•
possible	Seam	•	•
	Loose thread		

			Product for:		
Alteration			Thorax	Arm	
Pantyhose				_	
Panty area					
	Reduce circumference of pantyhose top				
	Incorporate U-part retrospectively				
	Shorten pantyhose top				
	Incorporate zip fastening	•			
	Enlarged open crotch				
	Incorporate stoma opening in pantyhose top				
	Waist edging				
	Gusset				
AT to AG wit	th garter belt	•			
AT to CT, ET,		•			
	th silicone top band	•			
_eg/arm	'				
Shorten sto	cking (chain-stitched edging)	•		•	
	band/waist band/strap	-		_	
	Sew in new sections of silicone top band				
	Sew on new silicone top band/waist band				
	Silicone top band up to 5 cm narrower or wider				
	Alter garter belt to AG with silicone top band	•		•	
	Alter shoulder support to CG with silicone top band				
	New strap for garter belt/ shoulder support (wear)				
Sew garter	belt onto leg stockings up to G	•			
Sew should	er support onto CG		_	•	
_ymph pad	/lining	•		•	
Finishing		•		•	
Foot (per fo	oot)				
Toe alteration	ons/foot length alteration	•			
Reduce cA	dimension	•			
Detach little	e toe from toe cap	•			
Hand altera	ations				
Shorten find	gers (thumb only)			•	
Shorten/tak	ke in glove			•	
	/alter lymph pad retrospectively			•	
Thorax	,			_	
	gthen, shorten, or narrower at the front/back)		•		
Arm sleeve	(per arm)		•		
_ower neck	line		•		
Shorten the	ırax		•		
Finishing			•		
	odifications				
Seam repai		•	•	•	
Run repair		•	•		
Hole repair					
	zip fastening				
· · · · · · · · · · · · · · · · · · ·					
Remove zip	-	•	•	•	
Repair zip f	the state of the s				


Please do **not** try on the stockings – redyeing is no longer possible after contact with skin
 Only possible with temporary seam

SIZING SYSTEMS

Sizing system V6+ VenoTrain impuls



Sizing system VenoTrain soft Arm



Length	short	long	Туре
€D	39 – 44	44,1 - 50	CG

ORDER PADS

We have developed a range of customer order pads for ordering standard and custom-made solutions. With these pads, the required version of a particular VenoTrain product can be ordered quickly and easily. We are happy to send you these order pads free of charge.

The different order sheets can also be downloaded at any time:

www.bauerfeind.com/downloads

STANDARD ORDE	Item number	
	VenoTrain	27120200000005
CUSTOM-MADE O	Item number	
Circular knit product	27120200001003	
Flat knit product	VenoTrain soft Arm	27112203021001
	VenoTrain delight leg	27120200000003
	VenoTrain delight arm	27120200001005
	VenoTrain delight glove	27120200001006
	VenoTrain delight toe caps	27120200001007
	VenoTrain curaflow leg	27142207001001
	VenoTrain curaflow toes	27142207001002
	VenoTrain curaflow arm/hand	27142207001003
	VenoTrain curaflow thorax	27141203001003

BAUERFEIND AG TERMS AND CONDITIONS OF SALE AND DELIVERY (GTC)

VERSION: FEBRUARY 2016

Prices are quoted as net prices in euros ex works and do not include statutory VAT.

Shipping is at the risk and on the account of the purchaser, even if the purchaser does not bear the costs of transportation. Disruptions of operations affecting us or any of our suppliers, lack of suitable personnel, strike, lockout, and force majeure are causes for extending delivery and relieve us of any delivery obligation for the duration of the obstruction.

Shipping costs are calculated as a flat rate of EUR 7 for a net total order value of up to EUR 50, and EUR 4 for a net total order value of up to EUR 200. No shipping costs are payable for custom-made solutions. Special shipping requests are calculated on a time and material basis.

Returns are only accepted outside the warranty terms if we have been notified in advance and approved the return (notification) with inclusion of the shipping documents and under the assumption that the returned goods are in faultless condition. We reserve the right to reject any returned consignment. Except where covered by warranty, we charge for processing returns. For returns within one month of the date of delivery, this charge is EUR 10.00 for each return consignment. For goods returned one month or more after the date of delivery, a handling fee of 20 percent of the net order value, or a maximum of EUR 25.00 per product and unit, is charged; for returns six to 12 months after the delivery date, the fee is 40 percent of the net order value, or a maximum fee of EUR 50.00 for each product and unit. Except where covered by warranty, no returns, including without limitation for custommade orders, will be accepted after 12 months following the date of delivery.

Payments are due immediately following delivery. We offer a four percent discount for payments made by direct debit and a three percent discount for payments made within 10 days of the invoice date. Otherwise payment is 30 days net. The buyer shall be considered in arrears once the 30 days from the invoice date have passed. From this point on, any further deliveries will only be made against advance payment. From the date that the customer is in arrears, we will charge default interest of eight percent above the applicable base interest rate. Payments by check, draft, etc. are deemed to be made only when the amount has been credited fully and irrevocably to our account.

Retention of proprietary rights

We retain the proprietary rights to all of the delivered goods until the purchaser has settled all of the outstanding amounts due to us as a result of the business relationship. The retention of proprietary rights shall also apply pro rata if the product in question has been processed or modified. If the goods are resold, all third-party claims arising therefrom up to the total amount owed to us shall be transferred to us without proviso. Notwithstanding the assignment, the buyer is entitled to the collection of these amounts. Our right to collect the amounts is

not affected by the preceding provision; however, we undertake to refrain from collection unless the buyer defaults on payment. If the buyer defaults on payment, we can demand that the buyer discloses the assigned claims and the names of the debtors, provides all information required for collection of the debt, delivers all required documents, and informs the debtors (third parties) of the assignment. The buyer shall not be entitled to pledge or assign the goods as security. Pledges, confiscation, or any other dispositions taken by a third party shall be promptly communicated to us by the buyer and all information and documents necessary for us to defend our rights shall be provided.

We undertake to release the security due to us at the seller's request to the extent that its value does not exceed that of the claim to be secured by more than 20 percent and provided the claim in question has not yet been settled.

Complaints

Warranty claims in respect of obvious defects shall be void if we do not receive immediate written notification of these claims in accordance with Section 377 of the German Commercial Code (HGB). In all other cases, the warranty period for warranty claims is one year after delivery of the purchased goods. If the goods are defective, we may replace or repair the goods at our discretion. For reasons of hygiene, repairs can only be made to goods that have been washed or which are in perfect sanitary condition. Upon receipt of unwashed goods, we reserve the right to claim a suitable sum for additional work. Provided that we fulfill our supplementary performance obligations, especially in regard to rectifying defects or delivering a defect-free product, the customer has no right to demand a price reduction or to withdraw from the contract, unless we fail to deliver our supplementary performance. Defects which are due to normal wear and tear, overuse, misuse, neglect of care, or noncompliance with the instructions for use are not covered by warranty.

Liability for damages is excluded unless there has been grossly negligent or willful misconduct on our part or on the part of our representatives or agents. This limitation of liability does not apply to damage due to injury of life, body, or health or the violation of a material contractual duty.

Place of jurisdiction

The place of performance for the delivery is the registered office of our company. The place of jurisdiction for all disputes arising from the business relation is Gera, Germany.

Severability clause

In the event that any individual provision of these terms and conditions of sale and delivery is or becomes invalid, the validity of the remaining provisions shall remain unaffected. The invalid provision shall be replaced with a valid provision which reflects the original economic intention of the parties as closely as possible.

BETTER QUALITY OF LIFE FOR YOUR CUSTOMERS!

Have you got any questions, suggestions, or requests? Please contact us:

BAUERFEIND AG

Triebeser Straße 16 07937 Zeulenroda-Triebes Germany

P +49 (0) 36628 66-20 00

F +49 (0) 36628 66-2999

E info@bauerfeind.com

AUSTRIA

Bauerfeind Ges.m.b.H. Hainburger Straße 33 1030 Wien

P +43 (0) 800 4430-130

F +43 (0) 800 4430-131

E info@bauerfeind.at

BENELUX

Bauerfeind Benelux B.V. Waarderveldweg 1 2031 BK Haarlem

THE NETHERLANDS

P +31 (0) 23 531-9427

F +31 (0) 23 532-1970

E info@bauerfeind.nl

BELGIUM

P +32 (0) 2 527-4060

F +32 (0) 2 792-5345

E info@bauerfeind.be

BOSNIA AND HERZEGOVINA

Bauerfeind d.o.o. Meše Selmovića 19 71000 Sarajevo

P +387 (0) 33 710-100

F +387(0)33619-422

E info@bauerfeind.ba

CROATIA

Bauerfeind d.o.o. Goleška 20 10020 Zagreb

P +385 (0) 1 6542-855

F +385 (0) 1 6542-860

E info@bauerfeind.hr

FRANCE

Bauerfeind France S.A.R.L. B.P. 50258 95957 Roissy CDG Cedex

P +33 (0) 1 4863-2896

F +33 (0) 1 4863-2963

E info@bauerfeind.fr

ITALY

Bauerfeind AG Filiale Italiana Piazza Don Enrico Mapelli 75 20099 Sesto San Giovanni (MI)

P +3902 8977 6310

F +390289775900

E info@bauerfeind.it

NORDIC

Bauerfeind Nordic AB Kungsholmstorg 8 11221 Stockholm

P +46 (0) 774 100-020

F +46 (0) 774 100-021

E info@bauerfeind.se

REPUBLIC OF MACEDONIA

Bauerfeind Dooel Skopje 50 Divizija 24 a 1000 Skopje

P +389(0)23179-002

F +389(0)23179-004

E info@bauerfeind.mk

SERBIA

Bauerfeind d.o.o. 102 Omladinskih brigada 11070 Novi Beograd

P +381 (0) 11 2287-050

F +381 (0) 11 2287-052

E info@bauerfeind.rs

SINGAPORE

Bauerfeind Singapore Pte Ltd. Blk 41 Cambridge Road #01-21 Singapore 210041

P +65 6396-3497

F +65 6295-5062

E info@bauerfeind.com.sq

SLOVENIA

Bauerfeind d.o.o. Dolenjska cesta 242 b 1000 Ljubljana

P +386 (0) 1 4272-941

F +386(0)1 4272-951

E info@bauerfeind.si

SPAIN

Bauerfeind Ibérica, S.A. Historiador Diago, 13 46007 Valencia

P +34 96 385-6633

F +34 96 385-6699

E info@bauerfeind.es

SWITZERLAND

Bauerfeind AG Vorderi Böde 5 5452 Oberrohrdorf

P +41 (0) 56 485-8242

F +41 (0) 56 485-8259

E info@bauerfeind.ch

UNITED ARAB EMIRATES

Bauerfeind Middle East FZ-LLC GF-09, District 1, 71 MPC Building Dubai Healthcare City P.O. Box 505116 Dubai

P +971 4 4335-685

F +971 4 4370-344

E info@bauerfeind.ae

UNITED KINGDOM

Bauerfeind UK 85 Tottenham Court Road London W1T 4TQ

P +44(0)121 446-5353

F +44(0)797457-5488

E info@bauerfeind.co.uk

USA

Bauerfeind USA, Inc. 1230 Peachtree St NE Suite 1000

Atlanta, GA 30309

P +1 800 423-3405

P +1 404 201-7800

F +1 404 201-7839

E info@bauerfeindusa.com